

# Lipase DC\* FS\*\*

Diagnostic reagent for quantitative in vitro determination of lipase in serum or plasma on photometric systems

## Order Information

Cat. No.	Kit size					
1 4321 99 10 021	R1	5 x	20 mL	+ R2	1 x	25 mL
1 4321 99 10 023	R1	1 x	800 mL	+ R2	1 x	200 mL
1 4321 99 10 930	R1	4 x	20 mL	+ R2	2 x	10 mL

## Summary

Lipases are enzymes which hydrolyze glycerol esters of long fatty acids. The enzyme and its cofactor colipase are produced in the pancreas, lipase being also secreted in small amounts by the salivary glands as well as by gastric, pulmonary and intestinal mucosa. Bile acids and colipase form micellar complexes with the lipids and bind lipase on the substrate/water interface. Determination of lipase is used for investigation of pancreatic disorders. In acute pancreatitis the lipase concentrations rise to 2 - 50 fold the upper reference limit within 4 – 8 hours after the beginning of abdominal pain peaking at 24 hours and decrease within 8 to 14 days. Elevated lipase values can also be observed in chronic pancreatitis and obstruction of the pancreatic duct. [1,2,3,4]

## Method

Enzymatic color test

A synthetically produced lipase substrate (1,2-o-dilauryl-rac-glycero-3-glutaric acid-(6-methylresorufin) ester) is added to a micro-emulsion which is specifically split by lipase in the presence of colipase and bile acids. The combination of lipase and bile acids make this specific and reliable for pancreatic lipase without any reaction due to lipolytic enzymes or esterases. The reagent composition has been thoroughly optimized so there are no serum matrix effects.

The generated methylresorufin-ester is spontaneously degraded to methylresorufin. The absorbance by this red dye is directly proportional to the lipase activity in the sample. [5,6,7]

## Principle

Lipase catalyses the reaction

1,2-o-Dilauryl-rac-glycero-3-glutaric acid(6-methylresorufin) ester

← Lipase / Colipase →

1,2-o-Dilauryl-rac-glycerin + Glutaric acid-(6-methylresorufin)-ester

Glutaric acid-(6-methylresorufin)-ester ← spontaneous degradation →

Glutaric acid + Methylresorufin

The increase in absorbance is determined photometrically.

## Reagents

### Components and Concentrations

#### Reagent 1:

Good's Buffer	pH 8.0	50 mmol/L
Taurodesoxycholate		4.3 mmol/L
Desoxycholate		8.0 mmol/L
Calcium chloride		15 mmol/L
Colipase (porcine)		2.2 mg/L

#### Reagent 2:

Tartrate Buffer	pH 4.0	7.5 mmol/L
Taurodesoxycholate		17.2 mmol/L
Color Substrate		≤ 0.65 mmol/L

### Storage Instructions and Reagent Stability

The reagents are stable up to the end of the indicated month of expiry, if stored at 2 – 8°C and contamination is avoided. Do not freeze the reagents and store them protected from light!

**Note:** A slight red precipitate may occur in reagent 2 which does not affect the performance of the test. Please do not resuspend before use!

### Waste Management

Please refer to local legal requirements.

### Reagent Preparation

The reagents are ready to use. Do not shake!

### Materials required but not provided

NaCl solution 9 g/L

General laboratory equipment

### Warnings and Precautions

1. Reagent 2: Warning. H319 Causes serious eye irritation. P280 Wear protective gloves/protective clothing/ eye protection/face protection. P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337+P313 If eye irritation persists: Get medical advice/attention.
2. Reagent 1 contains sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
3. Reagent 1: contains animal material. Handle the product as potentially infectious according to universal precautions and good clinical laboratory practices.
4. Many other clinical reagents contain lipase or high concentrations of detergents. Avoid contamination and carry over! Special care should be taken in combination with triglycerides, HDL and LDL reagents. Cuvettes and other glassware must be cleaned thoroughly after being used for other assays. In case of automated measurement refer to the instrument manual for special washing programs.
5. In very rare cases, samples of patients with gammopathy might give falsified results [11].
6. Please refer to the safety data sheets and take the necessary precautions for the use of laboratory reagents. For diagnostic purposes, the results should always be assessed with the patient's medical history, clinical examinations and other findings.
7. For professional use only!

## Specimen

Serum or heparin plasma

Stability [8]:	7 days	at	20 – 25°C
	7 days	at	4 – 8°C
	1 year	at	-20°C

Discard contaminated specimens! Only freeze once!

## Assay Procedure

**Application sheets for automated systems are available on request.**

Wavelength	580 nm, Hg 578 nm
Optical path	1 cm
Temperature	37°C
Measurement	Against reagent blank

	Blank	Sample
Sample or calibrator	-	20 µL
Dist. water	20 µL	-
Reagent 1	1000 µL	1000 µL
Mix carefully (do not shake), incubate 1 to 5 min. Start reaction by adding reagent 2:		
Reagent 2	250 µL	250 µL
Mix, incubate 2 min at 37°C, read absorbance and start stop watch. After exactly 1 and 2 min read absorbance again and then calculate ΔA/min.		

$$\Delta A/\text{min} = [\Delta A/\text{min sample or calibrator}] - [\Delta A/\text{min blank}]$$

## Calculation

With calibrator:

$$\text{Lipase [U/L]} = \frac{\Delta A / \text{min Sample}}{\Delta A / \text{min Calibrator}} \times \text{Conc. Calibrator [U/L]}$$

## Conversion factor

$$\text{Lipase [U/L]} \times 0.0167 = \text{Lipase [\mu kat/L]}$$

## Calibrators and Controls

DiaSys TruCal U calibrator is recommended for the calibration of automated photometric systems. The assigned values of the calibrator have been made traceable to the molar extinction coefficient of an available measuring method. DiaSys TruLab N and P controls should be assayed for internal quality control. Use of human based controls is strictly recommended. Each laboratory should establish corrective action in case of deviations in control recovery.

	Cat. No.	Kit size
TruCal U	5 9100 99 10 063	20 x 3 mL
	5 9100 99 10 064	6 x 3 mL
TruLab N	5 9000 99 10 062	20 x 5 mL
	5 9000 99 10 061	6 x 5 mL
TruLab P	5 9050 99 10 062	20 x 5 mL
	5 9050 99 10 061	6 x 5 mL

## Performance Characteristics

### Data evaluated on BioMajesty® JCA-BM6010/C

Exemplary data mentioned below may slightly differ in case of deviating measurement conditions.

Measuring range up to 300 U/L. When values exceed this range samples should be diluted 1 + 1 with NaCl solution (9 g/L) and the result multiplied by 2.	
Limit of detection***	5 U/L
Interfering substance	Interferences < 10% up to
Ascorbic acid	60 mg/dL
Bilirubin	70 mg/dL
Ditaurobilirubin	60 mg/dL
Hemoglobin	600 mg/dL
Lipemia (triglycerides)	2000 mg/dL
N-acetylcysteine (NAC)	2000 mg/L
For further information on interfering substances refer to Young DS [9].	

Precision			
Within run (n=20)	Sample 1	Sample 2	Sample 3
Mean [U/L]	30.9	60.9	286
CV [%]	1.26	0.611	0.263
Total Precision (n=80)	Sample 1	Sample 2	Sample 3
Mean [U/L]	30.2	59.9	284
CV [%]	2.01	1.20	1.10

Method comparison (n=107)	
Test x	Competitor Lipase (Cobas c311)
Test y	DiaSys Lipase DC FS (BioMajesty® JCA-BM6010/C)
Slope	0.982
Intercept	-0.168 U/L
Coefficient of correlation	0.999

\*\*\* according to CLSI document EP17-A2, Vol. 32, No. 8

## Reference Range [10]



$$\leq 60 \text{ U/L} \quad \leq 1.00 \mu \text{kat/L}$$

Each laboratory should check if reference ranges are transferable to its own patient population and determine own reference ranges if necessary.

## Literature

- Lorentz K. Lipase. In: Thomas L, editor. Clinical laboratory diagnostics. 1<sup>st</sup> ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 95-7.
- Moss DW, Henderson AR. Digestive enzymes of pancreatic origin. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3<sup>rd</sup> ed. Philadelphia: W.B Saunders Company; 1999. p. 689-708.
- Tietz N, Shuey DF. Lipase in serum – the elusive enzyme: an overview. Clin Chem 1993; 39: 746-56.
- Lott J, Patel ST, Sawhney AK, Kazmierczak SC, Love JE. Assays of serum lipase: analytical and clinical considerations. Clin Chem 1986; 32: 1290-1302.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv Clin Enzymol 1986; 4: 60-7.
- Borgström B. The action of bile salts and other detergents on pancreatic lipase and the interaction with colipase. Biochimica et Biophysica Acta 1977; 488: 381-91.
- Gargouri Y, Julien R, Bois A, Verger R, Sarda L. Studies on the detergent inhibition of pancreatic lipase activity. J of Lipid Research 1983; 24: 1336-42.
- Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1<sup>st</sup> ed. Darmstadt: GIT Verlag; 2001; p. 36-7.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
- Junge W, Abicht K, Goldman J. Evaluation of the colorimetric liquid assay for pancreatic lipase on Hitachi analyzers in 7 clinical centres in Europe. Clin Chem Lab Med 1999; 37, Special suppl: 469.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanism, detection and prevention. Clin Chem Lab Med 2007; 45(9): 1240-1243.

## Manufacturer

  DiaSys Diagnostic Systems GmbH  
Alte Strasse 9 65558 Holzheim Germany